


1079 INJECTOR W\CO2 CRYO

Subassembly

03-925445-11 (120v CO2)

03-925445-12 (230v CO2)

03-925445-13 (101v CO2)




TOOLS LIST

1. T-3734 ASSY. FIX.

OTHER PARTS NOT SHOWN.

03-925572-01 CLAMP, TUBING



OPERATION SEQUENCE

NOTE:

1. HEATER CARTRIDGE & PROBE CONTAINED BY CLAMSHELLS. END OF HEATER CARTRIDGE TO BE FLUSH WITH THE BOTTOM OF CLAMSHELLS.
2. CRYO TUBE & CLAMSHELL ASSY TO BE MOUNTED AS FAR UP THE CHAMBER TUBE AS POSSIBLE.