


1079 INJECTOR WITH LN2 CRYO

Subassembly

03-925445-21 (120v LN2)

03-925445-22 (230v LN2)

03-925445-23 (101v LN2)




TOOLS LIST

1. T-3734
ASSY. FIX.

OTHER PARTS NOT SHOWN.

03-925572-01
CLAMP, TUBING



OPERATION SEQUENCE

NOTE:

1. HEATER CARTRIDGE & PROBE CONTAINED BY CLAMSHELLS. END OF HEATER CARTRIDGE TO BE FLUSH WITH THE BOTTOM OF CLAMSHELLS.
2. CRYO TUBE & CLAMSHELL ASSY TO BE MOUNTED AS FAR UP THE CHAMBER TUBE AS POSSIBLE.