



1079/Manual Pneu. Installation



03-918532-00
VAR SEPTUM PURGE



12-901163-00
SCRW., SEMS 8-32 X 3/8 LG. EXT. STL 4EA

[CLICK HERE FOR
HEATER CONNECTION
TO J41 ON POWER BD.](#)

03-925425-01
120V HTR. CART.

03-925425-02
230V HTR. CART

03-925425-03
101V HTR. CART.

REF

03-925516-01
S/A 1079 INJ. MOD. MANUAL 120v

03-925516-02
S/A 1079 INJ. MOD. MANUAL 230v

03-925516-03
S/A 1079 INJ. MOD. MANUAL 101v


[CLICK HERE FOR
PROBE CONNECTION
TO J7 ON MAIN BOARD](#)

REFERENCE
03-925374-01
TEMPERATURE PROBE

PARTS NOT SHOWN

28-694581-00
FERRULE # 100 / 0.5

03-949551-00
NUT, CAPILLARY, INJ



NOTE: DETAIL OF SOLENOID VALVE

[CLICK HERE](#) TO SEE SOLENOID CONNECTION TO J42 ON POWER BD.

